# 

Latiful Kabir

## Contents

1	Resources	3
2	Quick start on basestar	4
3	Event length and file size calculation	5
4	The data file structure	6
5	The ADC channel to wire map	9
6	Setting up local version of n3He analysis library on basestar	10
7	Setting up local version of data browser on basestar	12
8	Current tree structure in n3He analysis libary	14
9	Sample analysis	16
10	Reference for TTreeRaw class	20

#### 1 Resources

All the startup software and manual related to  $n^3He$  experiment can be found from the official git repository with detailed instruction.

The easiest way is to go to n3He wiki (n3he.wikispaces.com) and then click on software from the left panel.

Alternatively, here is a direct link.

Any new change goes to this git repository. You might want to clone the entire repository and pull periodically to be updated. Or you can also download the last release of only what your are interested in.

### 2 Quick start on basestar

On basestar the data is being transferred and saved to the directories:

```
/mnt/idata01/data/ (run number:1-23661),
/mnt/idata02/data/ (run number:23662 - ),
/mnt/idata03/data/ (run number: - )
```

and the analysis library is compiled in a shared directory /home/npdg/n3He/libn3He/lib

So a quick start using the compiled library can be as follows from any user account:

1. Add the following lines to the .bashrc file & save it

- 2. Start a new terminal, go to /home/npdg/n3He/libn3He/analysis/ directory & try running sample analysis scripts from ROOT.
- 3. The data browser GUI (named as n3HeData) can be opened issuing the command/home/npdg/n3He/n3HeData/n3HeData from the terminal. Copy the binary to your home directory if you will be using the GUI frequently.

Other user specific customization can be achieved following the instruction in the ReadMe file in the respective directory.

For a local version of the library and data browser please read the corresponding section.

### 3 Event length and file size calculation

The event length set in the clean DAQ at 50 KHz sample rate is :830 . Theoretically the maximum possible value is  $50 \text{KHz} \times 16.66 \text{ ms} = 833 \text{ samples}$  per T0 .

But 833 event length gives occasional overlap. So we set to 830.

For dirty DAQ at 100 KHz the theoretical number of samples 100 KHz x 16.66 ms = 1666

But to avoid overlap we set to 1660.

More over the DAQ has fixed dead time(readout time) of 35 samples(with no averaging) at the end of any event. This amount of time will be missed for every event.

Clean DAQ event length 830 with nacc=16,16 with hi resolution mode=1 Dirty DAQ event length 1660 with nacc=1,1 with hi resolution mode=0 where nacc=n,n indicates how many samples being averaged.

Thus number of sample per event:

Clean DAQ:  $(830-35)/16=49.68 \sim 50 \text{ (1 header} + 49 \text{ samples)}$ 

Dirty DAQ: (1660-35)=1625 (1 header + 1624 samples)

Run Length calculation:

With 25000 T0 per run:

Clean DAQ file size: 25000 T0 x 50 samples x 4 Byte per sample x 48 Channels  $=240 \times 10^6$  Bytes

Dirty DAQ file size (before process): 25000 T0 x 1625 samples x 4 Bytes per sample x 8 channels =  $1300 \times 10^6$  Bytes

Dirty DAQ file seize (after process) : 25000 T0 x 1625 samples x 4 Bytes per sample x 2 channels =  $325 \times 10^6$  Bytes

#### 4 The data file structure

48 Clean DAQ channels divided into two modules:

Each sample is 4 bytes(in hexdump one contiguous pair consists one sample or 4 bytes). Out of this 32 bit(4 bytes), our data is 24 bit and remaining least significant 8 bits are channel ID.



Figure 1: Typical view of hexdump

The above hexdump to be interpreted as follows: With: mod= module ch = channel

```
      mod1event1sample1Ch0
      mod1event1sample1ch1
      mod1event1sample1ch8

      mod1event1sample1Ch9
      mod1event1sample1ch10
      mod1event1sample1ch11
      mod1event1sample1ch16

      mod1event1sample1Ch17
      mod1event1sample1ch18
      mod1event1sample1ch19
      mod1event1sample1ch23

      mod2event1sample1Ch0
      mod2event1sample1ch1
      mod2event1sample1ch8

      mod2event1sample1Ch9
      mod2event1sample1ch10
      mod2event1sample1ch16

      mod2event1sample1Ch17
      mod2event1sample1ch18
      mod2event1sample1ch19
```

mod1eve	nt2sample1Ch0	mod1event2sample1ch1	mod1event2sample1ch3	mod1event2sample1c	:h8
mod1eve	nt2sample1Ch9	mod1event2sample1ch1	.0 mod1event2sample1ch11	modlevent2sample	1ch16
mod1eve	nt2sample1Ch17	mod1event2sample1ch	118 mod1event2sample1ch19	mod1event2sampl	.e1ch23
mod2eve	nt2sample1Ch0	mod2event2sample1ch1	mod2event2sample1ch3	mod2event2sample1c	:h8
mod2eve	nt2sample1Ch9	mod2event2sample1ch1	.0 mod2event2sample1ch11	mod2event2sample	1ch16
mod2eve	nt2sample1Ch17	mod2event2sample1ch	118 mod2event2sample1ch19	mod2event2sampl	.e1ch23

Up to N number of events.

Now the first sample of any event is the event header with following structure:

mod1event1sample1Ch0=mod1event1sample1Ch1=mod1event1sample1Ch3 = Event Signature-1 (0xaa55f154),

mod1event1sample1Ch4= Event Number

mod1event1sample1Ch5 = checksum using path-1

mod1event1sample1Ch6 = sample number

mod1event1sample1Ch7 = checksum using path-2

Then this pattern repeats 3 more times (i.e. in quanta of 8 channels) up to channel-23

mod2event1sample1Ch0 = mod2event1sample1Ch1 = mod2event1sample1Ch3 = Event Signature-2 (0xaa55f15f)

mod1event1sample1Ch4= 0 (always)

mod2event1sample1Ch5 = checksum using path-1

mod2event1sample1Ch6 = sample number

mod2event1sample1Ch7 = checksum using path-2

Then this pattern repeats 3 more times (i.e. in quanta of 8 channels) up to channel-23

For Dirty DAQ the data is taken in 8 channels (bank mask B) with one module only and then processed to 2 channels. On Batch panel, M1 signal is connected to marked channel-26 and RFSF signal is connected to marked channel 27. This corresponds to ADC channel-5 (with checksum)

and channel-6 (with sample number) where for ADC channel number starts with 0.				
This for Dirty DAQ after the processing,				
event1sample1Ch0 event1sample1ch1				
1				
event2sample1Ch0 event2sample1ch1				
Up to N events.				
with the first sample of any event being checksum and sample number i.e.				
event1sample1ch0 = checksum.				
event1sample1ch1 = sample number.				
eventisampietent – sampie number.				

## 5 The ADC channel to wire map

Out of 48 clean ADC channels, only 36 ADC channels are connected to the pre-amps. Eventually the ADC channel to detector wire map to be handled by the library. However following is the map for reference:

```
Number of layers = 16;
Number of wires per layer= 9;
Layer_to_DAQ_map[Nlayers]={21, 23, 21, 23, 21, 23, 21, 23,
                       22, 24, 22, 24, 22, 24, 22, 24};
Layer_to_ADC_channel_map[16][9] =
     \{0,1,2,3,4,5,6,7,8\},\
     \{0,1,2,3,4,5,6,7,8\},
     {9,10,11,12,13,14,15,16,17},
     {9,10,11,12,13,14,15,16,17},
     {24,25,26,27,28,29,30,31,32},
     {24,25,26,27,28,29,30,31,32},
     {33,34,35,36,37,38,39,40,41},
     {33,34,35,36,37,38,39,40,41},
     \{0,1,2,3,4,5,6,7,8\},
     \{0,1,2,3,4,5,6,7,8\},
     {9,10,11,12,13,14,15,16,17},
     {9,10,11,12,13,14,15,16,17},
     {24,25,26,27,28,29,30,31,32},
     {24,25,26,27,28,29,30,31,32},
     {33,34,35,36,37,38,39,40,41},
     {33,34,35,36,37,38,39,40,41},
   };
```

## 6 Setting up local version of n3He analysis library on basestar

Eventually you will want to set up your own version of the analysis library and ROOT environment. This way you can modify any part of the library and add more functionality.

- 1. Download the source code from here. lib3He is the analysis library for n3He experiment.
- 2. Make any necessary changes in Constants.h file that is required.
- 3. Do make to compile the library.
- 4. This will produce libn3He.so (shared library will be inside lib directory).
- 5. Place the .so file in a directory under LD\_LIBRARY\_PATH .
- 6. Now start root and load the Library as: gSystem- >Load("libTree") & gSystem- >Load("libn3He.so"). (For Online analysis)
- 7. For analysis from a script if you include TTree.h file then you need not to do gSystem— >Load("libTree"); Just load gSystem— >Load("libn3He.so"). You need to give full path unless the directory is included in LD\_LIBRARY\_PATH.
- 8. If you put the rootlogon.C file in macros directory under Root installation directory, then the library will be loaded automatically and step-6 is NOT necessary.
- 9. Now from your root script create a Tree by calling: TTreeRaw \*my\_tree = new TreeRaw(runNumber#) or Just TTreeRaw t(runNumber#)
- 10. Do  $my\_tree > Print()$  to print the tree and branch structure.
- 11. Now do what ever analysis you want using my\_tree.
- 12. Try running example analysis scripts in "analysis" directory.
- 13. To make life easier it's convenient to put the following command into your  $\sim$ /.bash\_profile or  $\sim$ /.bashrc file:

Note: This version of the library works both for ROOT 5 and ROOT 6.

## 7 Setting up local version of data browser on basestar

- 1. Download the source code from here.
- 2. In bin directory: contains just binary files(obtained after doing make) named n3HeData.
- 3. In libn3He directory: Contains all the library required for running the Data check GUI.
- 4. Modify and compile the library (Unless you have already set up the library): You need to change the Data file directory from Constants.h in libn3He to appropriate directory. Before you do make do: make clean in the same directory. and then make a fresh shared binary files after you make any changes.
- 5. Place .so file under LD\_LIBRARY\_PATH: Now make sure the shared library (libn3He.so) file is in a directory under your LD\_LIBRARY\_PATH.

OR, a more professional way is as follows: Open /libn3He/bin/thisn3He.sh and define n3HeROOT path variable to the location where libn3He is located. for example: n3HeROOT=/home/siplu/libn3He/

Now include command in your .bashrc file to run thisn3He.sh file each time you open the terminal. i.e. include the following lines:

- 6. Produce binary for GUI: To produce a new binary file named n3HeData, go to n3HeData directory, open makefile and change LIB\_INCLUDE and GLIBS to appropriate location for you and then do make. It will produce n3HeData binary file in the same directory.
- 7. Run the GUI: Now run the binary file n3HeData by doing ./n3HeData from your recently compiled verson in /n3He\_DAQ\_GUI/n3HeData/ .

8. Modify the .desktop file (included one only for Ubuntu distribution) and .sh file in bin directory accordingly and place it in your desktop if you want to run the GUI just by double clicking from your desktop.




Figure 2: The n3He data browser

## 8 Current tree structure in n3He analysis libary

Currently n3He Tree has five branches corresponding to four clean DAQ and one dirty DAQ, following is the leaf list in the n3He Tree.

```
DAQ21_LEAF "h21[48]/I:d21[49][48]/I"

DAQ22_LEAF "h22[48]/I:d22[49][48]/I"

DAQ23_LEAF "h23[48]/I:d23[49][48]/I"

DAQ24_LEAF "h24[48]/I:d24[49][48]/I"

DAQ30_LEAF "h30[2]/I:d30[1624][2]/I"
```

where h used to indicate header and d used to indicate detector signal. 21 to 24 are clean DAQs and 30 is dirty DAQ. Clean DAQs contain detector signals. Dirty DAQ (after data processing) ADC channel-0 is monitor-1 signal and ACD channel-1 is RFSF signal.

```
Terminal - kabir@basestar:
File Edit View Terminal Go Help
[kabir@basestar ~]$ root -l
root [0] gSystem->Load("libTree")
(int)0
root [1] gSystem->Load("libn3He.so")
(int)0
root [2] TTreeRaw t(19005)
Reading clean dag file: /mnt/idata01/data/run-19005data-21
Reading clean daq file: /mnt/idata01/data/run-19005data-22
Reading clean daq file: /mnt/idata01/data/run-19005data-23
Reading clean daq file: /mnt/idata01/data/run-19005data-24
Reading dirty daq file: /mnt/idata01/data/run-19005data-30
root [3] t.GetListOfLeaves().Print()
Collection name='TObjArray', class='TObjArray', size=16
OBJ: TLeafI h21 h21[48]
OBJ: TLeafI
                 d21
                          d21[49][48]
OBJ: TLeafI
                 h22
                          h22[48]
OBJ: TLeafI
                 d22
                          d22[49][48]
                          h23[48]
d23[49][48]
h24[48]
 OBJ: TLeafI
                 h23
 OBJ: TLeafI
                 d23
 OBJ: TLeafI
                 h24
OBJ: TLeafI
                          d24[49][48]
                 d24
                          h30[2]
OBJ: TLeafI
                 h30
OBJ: TLeafI
                 d30
                          d30[1624][2]
root [4]
```

Figure 3: The n3He leaf list

Now the library always skips the first four or five events since those might NOT be reliable. As a result the number of events in any branch for a typical n3He run is 24996 or 24995(This offset is set dynamically).




Figure 4: The n3He tree structure

## 9 Sample analysis

Sample online analysis on basestar:

```
//OnlineAnalysis.C
//Demo Online Analysis using n3He Library.(By Online I mean 'from
    CINT, doing analysis on the fly, less thoughtful but preferred
    by in some conditions')
//Author: Latiful Kabir
//Date: 12/23/14

void OnlineAnalysis()
{
    gSystem->Load("libTree"); //You need to load libTree first in
        order to Load libn3He. This is not necessary if you include
        TTree.h
    gSystem->Load("libn3He.so");

    TTreeRaw *t=new TTreeRaw(19900);
    t->Draw("d21[][0]:Iteration$");
}
```

This script when you run using root -l OnlineAnalysis.C will produce the following output:



Figure 5: The Output from OnlineAnalysis.C

Sample offline analysis on basestar:

```
//OfflineAnalysis.C
//Demo Offline Analysis using n3He Library.(By Offline I mean 'in
    a script more thoughtful and serious analysis unlike from CINT)
//This script shows how to accress Tree using SetAddress
// and plots only the all event/pulses of channel-0
//Author: Latiful Kabir
//Date: 01/14/15
//This is the fastest and most preferred method for reading Tree
#include<TTree.h>
#include<TBranch.h>
#include<TGraph.h>
```

```
void OfflineAnalysis(){
   //Load the library unless loaded automatically by ROOT
    gSystem->Load("libTree");
    gSystem->Load("libn3He.so");
 //Create a TTreeRaw object with desired run number
 TTreeRaw *t=new TTreeRaw(17900);
 t->Print(); // Print to see what's inside the Tree
 int ch=0; //Channel to analyze
 //Create a struc buffer to keep your events
 struct myData
     int header[48];
     int det[49][48];
 };
 myData md;
 //Get the branch you want to analyze
 TBranch *b=t->GetBranch("b21");
 b->SetAddress(&md.header[ch]);
 TGraph *g=new TGraph();
 //Loop through all the events in the run.
 for(int i = 0;i < b->GetEntries();i++)
     //Load the samples for a event/pulse in buffer
     b->GetEntry(i);
     //Loops through the sample for the loaded event
     for(int k=0;k<49;k++)</pre>
    g->SetPoint(i*49+k,i*49+k,md.det[k][ch]);
 }
 g->Draw("AP");
```

```
delete t;
}
```

This script when you run using root -l OfflineAnalysis.C will produce the following output when zoomed in:




Figure 6: The output from OfflineAnalysis.C

## 10 Reference for TTreeRaw class

The base Class is TTree.

As a result all data and member functions from TTree are automatically inherited. For a list of TTree data and member functions go here

This list will be updated gradually as new functionality is added to the TTreeRaw class.